

For Research Use Only

CoraLite® Plus 647 Anti-Rat CD43 (W3/13) Mouse IgG2a Recombinant Antibody



Catalog Number: **CL647-65607**

Basic Information

Catalog Number:

CL647-65607

Size:

100ug, 500 µg/ml

Source:

Mouse

Isotype:

IgG2a

GenBank Accession Number:

GeneID (NCBI):

24796

Full Name:

sialophorin

Calculated MW:

38 kDa

Purification Method:

Protein A purification

CloneNo.:

W3/13

Excitation/Emission maxima
wavelengths:

654 nm / 674 nm

Applications

Tested Applications:

FC

Species Specificity:

rat

Background Information

Storage

Storage:

Store at 2-8°C. Avoid exposure to light.

Storage Buffer:

PBS with 0.09% sodium azide.

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1x10⁶ rat splenocytes were surface stained with 0.25ug CoraLite® Plus 647 Anti-Rat CD43 (W3/13) Mouse IgG2a Recombinant Antibody (CL647-65607, Clone: W3/13), or CoraLite® Plus 647 Mouse IgG2a Isotype Control (C1.18.4) (CL647-65208, Clone: C1.18.4). Cells were not fixed.



1x10⁶ rat splenocytes were surface stained with 0.25ug CoraLite® Plus 647 Anti-Rat CD43 (W3/13) Mouse IgG2a Recombinant Antibody (CL647-65607, Clone: W3/13)(red), or CoraLite® Plus 647 Mouse IgG2a Isotype Control (C1.18.4) (CL647-65208, Clone: C1.18.4) (blue). Cells were not fixed.