

For Research Use Only

CoraLite® Plus 647-conjugated NeuN Monoclonal antibody



Catalog Number:CL647-66836

Basic Information

Catalog Number:

CL647-66836

Size:

1 mg/mL

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG28016

GenBank Accession Number:

NM_001082575

GeneID (NCBI):

146713

Full Name:

hexaribonucleotide binding protein 3

Purification Method:

Protein G purification

CloneNo.:

3A4C1

Excitation/Emission maxima

wavelengths:

654 nm / 674 nm

Applications

Tested Applications:

FC (Intra)

Species Specificity:

Human, mouse, rat

Background Information

NeuN, encoded by FOX3, is a neuron-specific nuclear protein. Anti-NeuN stains exclusively neuronal cells in the central and peripheral nervous systems, especially postmitotic and differentiating neurons, as well as terminally differentiated neurons. Anti-NeuN has been used widely as a reliable tool to detect most postmitotic neuronal cell types. The immunohistochemical staining is primarily localized in the nucleus of the neurons with lighter staining in the cytoplasm.

Storage

Storage:

Store at -20°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1X10⁶ SH-SY5Y cells were intracellularly stained with 0.2 ug CoraLite® Plus 647 Anti-Human NeuN (CL647-66836, Clone:3A4C1) (red), or 0.2 ug APC-65128 ; APC Mouse IgG2b Isotype Control (APC-65128, Clone: MPC-11) (blue). Cells were fixed and permeabilized with Transcription Factor Staining Buffer Kit (PF00011).