

For Research Use Only

CoraLite® Plus 647-conjugated PCSK9 Recombinant antibody

Catalog Number: CL647-83004-4



Basic Information

Catalog Number:

CL647-83004-4

Size:

1000 µg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_153565.2

GeneID (NCBI):

100102

UNIPROT ID:

Q80W65

Full Name:

proprotein convertase
subtilisin/kexin type 9

Calculated MW:

75KD

Purification Method:

Protein A purification

CloneNo.:

230536E4

Recommended Dilutions:

IF/ICC 1:50-1:500

Excitation/Emission maxima
wavelengths:

654 nm / 674 nm

Applications

Tested Applications:

IF/ICC

Species Specificity:

mouse

Positive Controls:

IF/ICC : HepG2 cells,

Background Information

Storage

Storage:

Store at -20°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (-20°C Ethanol) fixed HepG2 cells using CoraLite® Plus 647 PCSK9 antibody (CL647-83004-4, Clone: 230536E4) at dilution of 1:100.