

For Research Use Only

CoraLite® Plus 750 Anti-Mouse CD90.2 (30-H12)



Catalog Number: CL750-65088

Basic Information

| | | |
|---------------------------------------|---|---|
| Catalog Number: CL750-65088 | GenBank Accession Number: BC054436 | Purification Method: Affinity purification |
| Size: 100ug, 0.5 mg/ml | GeneID (NCBI): 21838 | CloneNo.: 30-H12 |
| Source: Rat | Full Name: thymus cell antigen 1, theta | Excitation/Emission maxima wavelengths: 755 nm / 780 nm |
| Isotype: IgG2b, kappa | | |

Applications

Tested Applications:
FC

Species Specificity:
Mouse

Background Information

CD90 (Thy-1) is a 25 kDa, GPI-linked membrane glycoprotein that belongs to immunoglobulin superfamily (PMID: 6177036; 6153212). Originally described as a brain thymus cross-reactive antigen, it is found in large quantities on mouse and rat thymocytes and central nervous system cells (PMID: 83175). CD90 has been postulated to be involved in cellular recognition, adherence, and T cell activation (PMID: 7683034). Mouse CD90 (Thy-1) is defined as a differentiation alloantigen, represented as two serologically distinguishable allelic forms, designated CD90.1 (Thy-1.1) and CD90.2 (Thy-1.2) (PMID: 83175). CD90.2 is expressed by thymocytes, T cells, neurons and hematopoietic stem cells in most strains of mice (C3H, BALB/c, C57BL/6, DBA, and others), whereas CD90.1 is expressed in some mouse strains like AKR (PMID: 5919593).

Storage

Storage:
Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:
PBS with 0.09% sodium azide.

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1X10⁶ mouse thymocytes were surface stained with 0.5 ug CoraLite® Plus 750 Anti-Mouse CD90.2 (CL750-65088, Clone: 30-H12) (red), or 0.5 ug Control Antibody. Cells were not fixed.