

For Research Use Only

# CoraLite® Plus 750 Anti-Mouse CD41 (MWReg30)



Catalog Number: CL750-65189

## Basic Information

Catalog Number:

CL750-65189

Size:

100ug, 0.5 mg/ml

Source:

Rat

Isotype:

IgG1, kappa

GenBank Accession Number:

BC120493

GeneID (NCBI):

16399

UNIPROT ID:

Q9QUM0

Full Name:

integrin alpha 2b

Purification Method:

Affinity purification

CloneNo.:

MWReg30

Excitation/Emission maxima  
wavelengths:

755 nm / 780 nm

## Applications

Tested Applications:

FC

Species Specificity:

Mouse

## Background Information

### Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide.

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



1X10<sup>6</sup> Balb/c mouse peripheral blood platelets were surface stained with 0.2 ug CoraLite® Plus 750 Anti-Mouse CD41 (CL750-65189, Clone: MWReg30) (red), or 0.2 ug Isotype Control. Cells were not fixed.



1X10<sup>6</sup> Balb/c mouse peripheral blood platelets were surface stained with 0.2 ug CoraLite® Plus 750 Anti-Mouse CD41 (CL750-65189, Clone: MWReg30), or 0.2 ug Isotype Control. Cells were not fixed.