

For Research Use Only

# CoraLite® Plus 750 Anti-Mouse CD94 (18D3)



Catalog Number: CL750-65240

## Basic Information

Catalog Number:

CL750-65240

Size:

100ug, 0.5 mg/ml

Source:

Rat

Isotype:

IgG2a, kappa

GenBank Accession Number:

BC120854

GeneID (NCBI):

16643

UNIPROT ID:

O54707

Full Name:

killer cell lectin-like receptor,  
subfamily D, member 1

Purification Method:

Affinity purification

CloneNo.:

18D3

Excitation/Emission maxima  
wavelengths:

755 nm / 780 nm

## Applications

Tested Applications:

FC

Species Specificity:

Mouse

## Background Information

### Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide.

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



1x10<sup>6</sup> mouse splenocytes were surface stained with FITC Anti-Mouse CD49b (HMA2) (FITC-65245, Clone: HMA2) and 0.5 ug CoraLite® Plus 750 Anti-Mouse CD94 (CL750-65240, Clone: 18D3) or CoraLite® Plus 750 Rat IgG2a isotype control. Cells were not fixed.