

For Research Use Only

CoraLite® Plus 750-conjugated RILP Recombinant antibody

Catalog Number: CL750-82996-4



Basic Information

Catalog Number: CL750-82996-4	GenBank Accession Number: BC031621	Purification Method: Protein A purification
Size: 1000 µg/ml	GeneID (NCBI): 83547	CloneNo.: 230208G2
Source: Rabbit	UNIPROT ID: Q96NA2	Recommended Dilutions: WB 1:500-1:1000
Isotype: IgG	Full Name: Rab interacting lysosomal protein	Excitation/Emission maxima wavelengths: 755 nm / 780 nm
Immunogen Catalog Number: AG4472	Calculated MW: 401 aa, 41 kDa	
	Observed MW: 50 kDa	

Applications

Tested Applications: WB	Positive Controls: WB : HeLa cells, HepG2 cells
Species Specificity: human	

Background Information

Storage

Storage:
Store at -20°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:
PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with CL750-82996-4 (RILP antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.