

For Research Use Only

CoraLite® Plus 750-conjugated MTA3 Recombinant antibody

Catalog Number:CL750-83042-3



Basic Information

Catalog Number:

CL750-83042-3

Size:

1000 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG6400

GenBank Accession Number:

BC053631

GeneID (NCBI):

57504

UNIPROT ID:

Q9BTC8

Full Name:

metastasis associated 1 family,
member 3

Calculated MW:

68 kDa

Observed MW:

62 kDa, 68 kDa

Purification Method:

Protein A purification

CloneNo.:

230350D12

Recommended Dilutions:

WB 1:1000-1:8000

**Excitation/Emission maxima
wavelengths:**

755 nm / 780 nm

Applications

Tested Applications:

WB, FC (Intra)

Species Specificity:

human, mouse

Positive Controls:

WB : MCF-7 cells, Raji cells

Background Information

The metastasis-associated (MTA) family contains three members, MTA1, MTA2 and MTA3 (PMID:12167865). Metastasis-associated protein 3 (MTA3) is a cell-type specific component of the Mi-2-NURD transcriptional co-repressor complex. It has a role in maintenance of the normal epithelial architecture through the repression of SNAI1 transcription in a histone deacetylase-dependent manner, and thus the regulation of E-cadherin levels.

Storage

Storage:

Store at -20°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with CL750-83042-3 (MTA3 antibody) at dilution of 1:4000 incubated at room temperature for 1.5 hours.



1×10^6 U-251 cells were intracellularly stained with 0.8 μ g Coralite® Plus 750 Mta3 Recombinant Antibody (CL750-83042-3, Clone:230350D12)(red), or 0.8 μ g Isotype Control (blue). Cells were fixed and permeabilized with Transcription Factor Staining Buffer Kit (PF00011).