

Catalog Number: CM00228

产品信息

Catalog Number:
CM00228

CAS号:
3483-12-3


分子式:
 $C_4H_{10}O_2S_2$

主要靶点:
Others

主要通路:
其他

分子量:
154.25

溶解度:
H₂O:200 mg/mL (1296.60 mM), DMSO:125 mg/mL (810.37 mM), Sonification is recommended.



描述

DTT (DL-Dithiothreitol) is a reagent for maintaining SH groups in the reduced state.

储存

Powder: -20°C for 3 years | In solvent: -80°C for 2 years