

For Research Use Only

Thalidomide 4-fluoride

Catalog Number: CM01674

产品信息

Catalog Number:
CM01674

CAS号:
835616-60-9


分子式:
 $C_{13}H_9FN_2O_4$

主要靶点:
IRAK|Ligand for E3 Ligase

主要通路:
NF- κ B信号通路|免疫与炎症|PROTAC

分子量:
276.22

溶解度:
DMSO:199 mg/mL (720.44 mM)



描述

Thalidomide 4-fluoride (Cereblon ligand 4) is the Thalidomide-based Cereblon ligand. Thalidomide 4-fluoride can be used in the recruitment of CRBN protein. Thalidomide 4-fluoride (Cereblon ligand 4) can be connected to the ligand for IRAK4 protein by a linker to form PROTAC IRAK4 degrader-1.

储存

Powder: -20°C for 3 years | In solvent: -80°C for 2 years