

Catalog Number: CM01970

### 产品信息

**Catalog Number:**  
CM01970

**CAS号:**  
147127-19-3


**分子式:**  
C<sub>9</sub>H<sub>14</sub>N<sub>5</sub>O<sub>4</sub>P

**主要靶点:**  
Others

**主要通路:**  
其他

**分子量:**  
287.21

**溶解度:**  
DMSO:5 mg/mL (17.41  
mM),ultrasonic and warming and  
heat to 60°C



### 描述

(S)-Tenofovir is the less active S-enantiomer of Tenofovir which is a nucleotide reverse transcriptase inhibitor.

### 储存

Powder: -20°C for 3 years | In solvent: -80°C for 2 years