

For Research Use Only

# Amino-PEG3-C2-acid



www.ptgcn.com

Catalog Number: CM02064

## 产品信息

**Catalog Number:**  
CM02064

**CAS号:**  
784105-33-5

**分子式:**  
C<sub>9</sub>H<sub>19</sub>NO<sub>5</sub>


**主要靶点:**  
ADC Linker

**主要通路:**  
ADC

**分子量:**  
221.2509

**MDL NO:**  
MFCD21607085

**Pubchem ID:**  
11608245



## 描述

Amino-PEG3-C2-acid is a cleavable PEG ADC linker used in the synthesis of antibody-drug conjugates (ADCs).

## 储存

**储存条件:**

粉末	-20°C	3年
----	-------	----

液体	-80°C	1年
----	-------	----

**运输条件:**

Shipped in cold pack

For technical support and original validation data for this product please contact

T: 027-87531629

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## 参考文献

- 1 Kevin HAMBLETT, et al. Anti-her2 bipolaropic antibody-drug conjugates and methods of use. WO2019173911A1.
-