

For Research Use Only

BODIPY 493/503



www.ptgcn.com

Catalog Number: CM02294

产品信息

Catalog Number:
CM02294

CAS号:
121207-31-6


分子式:
 $C_{14}H_{17}BF_2N_2$

主要靶点:
Others

主要通路:
其他

分子量:
262.11

溶解度:
DMF:Soluble,Chloroform:Soluble,Me1
mg/mL(908.38 uM),DMSO:2
mg/mL(7.63 mM),Sonification is
recommended.



描述

BODIPY 493/503 is a lipophilic fluorescent probe with Ex/Em of 493/503 nm. BODIPY 493/503 localizes to polar lipids and can be used to label cellular neutral lipid contents and for live and fixed cell applications.

储存

Powder: -20°C for 3 years | In solvent: -80°C for 2 years

For technical support and original validation data for this product please contact

T: 027-87531629

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.