

For Research Use Only

MC-Val-Ala-PAB-PNP



www.ptgcn.com

Catalog Number: CM02508

产品信息

Catalog Number:
CM02508

CAS号:
1639939-40-4


分子式:
C₃₂H₃₇N₅O₁₀

主要靶点:
Others

主要通路:
其他

分子量:
651.66

溶解度:
DMSO:50 mg/mL (76.73 mM), Sonification is recommended.



体外活性

ADCs consist of an antibody attached to an ADC cytotoxin through an ADC linker.

描述

MC-Val-Ala-PAB-PNP is a cleavable ADC linker. MC-Val-Ala-PAB-PNP can be used in the synthesis of antibody-drug conjugates (ADCs).

储存

Powder: -20°C for 3 years | In solvent: -80°C for 2 years

For technical support and original validation data for this product please contact

T: 027-87531629

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.