

Catalog Number: CM02675

产品信息

Catalog Number:
CM02675

CAS号:
314728-85-3


分子式:
 $C_{14}H_{18}N_2O_2$

主要靶点:
iGluR

主要通路:
神经科学|离子通道

分子量:
246.3

溶解度:
DMSO:2.6 mg/mL (10.56 mM), H₂O:50 mg/mL (203.00 mM)



描述

Sunifiram is a potent nootropic agent.

储存

Powder: -20°C for 3 years | In solvent: -80°C for 2 years