

Catalog Number: CM03060

### 产品信息

**Catalog Number:**  
CM03060

**CAS号:**  
304481-60-5

**分子式:**  
C<sub>14</sub>H<sub>12</sub>N<sub>2</sub>O<sub>4</sub>


**主要靶点:**  
ATPase

**主要通路:**  
离子通道

**分子量:**  
272.26

**溶解度:**

DMSO:50mg/ml(183.65mM),Sonification is recommended.,Etanol:3mg/ml(11.02mM),Sonification is recommended.



### 描述

KM 91104 is a cell permeable inhibitor of V-ATPase. KM91104 specifically targets the interaction between V-ATPase subunit A3 and subunit B2.

### 储存

Powder: -20°C for 3 years | In solvent: -80°C for 2 years