

For Research Use Only

1-Methyl-7-nitroisatoic anhydride



www.ptgcn.com

Catalog Number: CM03148

产品信息

Catalog Number:
CM03148

CAS号:
73043-80-8


分子式:
 $C_9H_6N_2O_5$

主要靶点:
Others

主要通路:
其他

分子量:
222.15

溶解度:
DMSO:50 mg/mL (225.07 mM), Sonification is recommended.



描述

1-Methyl-7-nitroisatoic anhydride (1M7) is a reagent. It detects local nucleotide flexibility, for probing 2'-hydroxyl reactivity.

储存

Powder: -20°C for 3 years | In solvent: -80°C for 2 years

For technical support and original validation data for this product please contact

T: 027-87531629

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.