

Catalog Number: CM03625

产品信息

Catalog Number:
CM03625

CAS号:
2222768-84-3


分子式:
C₂₂H₂₃F₂N₃O₂S

主要靶点:
Hedgehog/Smoothed

主要通路:
干细胞|G蛋白偶联受体

分子量:
431.5

溶解度:
DMSO:80mg/mL
(185.40mM), Sonification is
recommended.



描述

Dynarrestin is a cytoplasmic dynein 1 and 2 aminothiazole inhibitor that reversibly inhibits dynein 1-driven microtubule sliding in vitro and intracellular dynein 1 and 2-dependent processes without affecting ATP hydrolysis and interfering with ciliogenesis.

储存

Powder: -20°C for 3 years | In solvent: -80°C for 2 years