

For Research Use Only

Carbidopa monohydrate



www.ptgcn.com

Catalog Number: CM03634

产品信息

Catalog Number:
CM03634

CAS号:
38821-49-7


分子式:
C₁₀H₁₆N₂O₅

主要靶点:
Decarboxylase|Aryl Hydrocarbon Receptor

主要通路:
其他|免疫与炎症|代谢

分子量:
244.24

溶解度:
DMSO:5.5 mg/mL (22.52 mM)



描述

Carbidopa monohydrate (S(-)-Carbidopa) is an inhibitor of DOPA decarboxylase, preventing the conversion of levodopa to dopamine.

储存

Powder: -20°C for 3 years | In solvent: -80°C for 2 years

For technical support and original validation data for this product please contact

T: 027-87531629

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.