

For Research Use Only

M1001



www.ptgcn.com

Catalog Number: CM04461

产品信息

Catalog Number:
CM04461

CAS号:
874590-32-6


分子式:
C₁₇H₁₇N₃O₂S

主要靶点:
HIF|HIF/HIF Prolyl-Hydroxylase

主要通路:
表观遗传|代谢|血管生成

分子量:
327.4

溶解度:
DMSO:250 mg/mL (763.59 mM)



靶点活性

HIF-2 α :667 nM(Kd)

描述

M1001 is a HIF-2 α agonist.

储存

Powder: -20°C for 3 years | In solvent: -80°C for 2 years

For technical support and original validation data for this product please contact

T: 027-87531629

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.