

For Research Use Only

Triprolidine hydrochloride monohydrate



www.ptgcn.com

Catalog Number: CM04509

产品信息

Catalog Number:
CM04509

CAS号:
6138-79-0


分子式:
C₁₉H₂₅ClN₂O

主要靶点:
Histamine Receptor

主要通路:
G蛋白偶联受体|免疫与炎症|神经科学

分子量:
332.87

溶解度:
DMSO:186.3 mM



描述

Triprolidine hydrochloride monohydrate is a histamine H1 antagonist used in allergic rhinitis, asthma.

储存

Powder: -20°C for 3 years | In solvent: -80°C for 2 years

For technical support and original validation data for this product please contact

T: 027-87531629

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.