

Catalog Number: CM04747

产品信息

Catalog Number:
CM04747

CAS号:
206645-99-0


分子式:
 $C_{66}H_{99}N_{23}O_{17}S_2$

主要靶点:
Complement System

主要通路:
免疫与炎症

分子量:
1550.77

溶解度:
DMSO:300 mg/mL (193.45 mM)



描述

Compstatin, a cyclic tridecapeptide, which was originally discovered from phage-display libraries, is a highly potent and selective C3 inhibitor.

储存

Powder: -20°C for 3 years | In solvent: -80°C for 2 years