

For Research Use Only

Fmoc-Asn(Trt)-Ser(ψ(Me,Me)pro)-OH



www.ptgcn.com

Catalog Number: CM04825

产品信息

Catalog Number:
CM04825

CAS号:
920519-33-1


分子式:
C₄₄H₄₁N₃O₇

主要靶点:
Intermediate and Others

主要通路:
Biochemical Reagent

分子量:
723.8123

MDL NO:
MFCD10001368



描述

Fmoc-Asn(Trt)-Ser(ψ(Me,Me)pro)-OH is used for solid phase synthesis of peptides.

储存

储存条件:

粉末	-20°C	3年
----	-------	----

液体	-80°C	1年
----	-------	----

运输条件:

Shipped in cold pack

For technical support and original validation data for this product please contact

T: 027-87531629

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.