

For Research Use Only

# Fmoc-Asp-NH<sub>2</sub>

Catalog Number: CM04826

## 产品信息

Catalog Number:  
CM04826


分子量:  
354.36

CAS号:  
200335-40-6

分子式:  
C<sub>19</sub>H<sub>18</sub>N<sub>2</sub>O<sub>5</sub>

主要靶点:  
Others

主要通路:  
其他



## 体外活性

ADCs consist of an antibody attached to an ADC cytotoxin through an ADC linker.

## 描述

Fmoc-Asp-NH<sub>2</sub> is a cleavable ADC linker. Fmoc-Asp-NH<sub>2</sub> can be used in the synthesis of antibody-drug conjugates (ADCs).

## 储存

Powder: -20°C for 3 years | In solvent: -80°C for 2 years

For technical support and original validation data for this product please contact

T: 027-87531629

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.