

Catalog Number: CM05466

产品信息

Catalog Number:
CM05466

CAS号:
2479-49-4


分子式:
 $C_{17}H_{10}O_9$

主要靶点:
Others

主要通路:
其他

分子量:
358.26

溶解度:
DMSO:45mg/mL
(125.6mM), Sonification is
recommended., H₂O:8 mg/mL
(22.33 mM), Sonification is
recommended.



描述

Benzophenonetetracarboxylic acid can be used in the preparation of high-performance polyimides and as a curing agent for epoxy resins.

储存

Powder: -20°C for 3 years | In solvent: -80°C for 2 years