

For Research Use Only

# MCU-i4



www.ptgcn.com

Catalog Number: CM05528

## 产品信息

**Catalog Number:**  
CM05528

**CAS号:**  
371924-24-2


**分子式:**  
C<sub>23</sub>H<sub>27</sub>N<sub>3</sub>O<sub>2</sub>

**主要靶点:**  
Others|Mitochondrial  
Metabolism

**主要通路:**  
其他|代谢

**分子量:**  
377.48

**溶解度:**  
DMSO:5 mM



## 描述

MCU-i4 is a novel negative modulator of the MCU, binding MICU1 and impairing muscle cell growth.

## 储存

Powder: -20°C for 3 years | In solvent: -80°C for 2 years

For technical support and original validation data for this product please contact

T: 027-87531629

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.