

Catalog Number: CM05554

产品信息

Catalog Number:
CM05554

CAS号:
478336-92-4


分子式:
C₁₂H₁₃NO₄

主要靶点:
Others

主要通路:
其他

分子量:
235.2

溶解度:
DMSO:23.5 mg/mL(100 mM)



靶点活性

MIF:7 μM

描述

ISO-1 (MIF Antagonist) is an inhibitor of the dopachrome tautomerase activity, blocks the activation of NF-κB and TNF-α secretion from LPS-treated macrophages.

储存

Powder: -20°C for 3 years | In solvent: -80°C for 2 years