

Catalog Number: CM05761

产品信息

Catalog Number:
CM05761

CAS号:
693241-54-2


分子式:
 $C_{22}H_{18}N_2O_5$

主要靶点:
Others

主要通路:
其他

分子量:
390.39

溶解度:
DMSO:70mg/ml(179.31mM), Sonification
is recommended.



描述

EG1, a novel specific inhibitor of Pax2 transcription activation, targets the DNA binding domain and inhibits embryonic kidney development.

储存

Powder: -20°C for 3 years | In solvent: -80°C for 2 years