

Catalog Number: CM05763

### 产品信息

**Catalog Number:**  
CM05763

**CAS号:**  
1632250-49-7


**分子式:**  
C<sub>22</sub>H<sub>20</sub>ClN<sub>7</sub>O

**主要靶点:**  
Others

**主要通路:**  
其他

**分子量:**  
433.9

**溶解度:**  
DMSO:4.34 mg/mL (10 mM), Sonification is recommended.



### 体内活性

PF-06446846 reduces plasma PCSK9 and total cholesterol levels in rats following oral dosing. PF-06446846 is highly selective for the inhibition of PCSK9 translation. The mechanism of action employed by PF-06446846 reveals a previously unexpected tunability of the human ribosome that allows small molecules to specifically block translation of individual transcripts.

### 描述

PF 06446846 is a potent and selective PCSK9 inhibitor

### 储存

Powder: -20°C for 3 years | In solvent: -80°C for 2 years