

For Research Use Only

SKF-34288 hydrochloride

Catalog Number: CM05831

产品信息

Catalog Number:
CM05831

CAS号:
320386-54-7


分子式:
 $C_6H_6ClNO_2S$

主要靶点:
Others

主要通路:
其他

分子量:
191.63

溶解度:
H₂O: 4 mg/mL (20.87 mM), Sonification is recommended.



描述

SKF-34288 hydrochloride is an inhibitor of phosphoenolpyruvate carboxykinase (PEPCK), and inhibits Asn metabolism and results in an increase in amino acids and amides.

储存

Powder: -20°C for 3 years | In solvent: -80°C for 2 years

For technical support and original validation data for this product please contact

T: 027-87531629

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.