

For Research Use Only

# Acalisib



www.ptgcn.com

Catalog Number: CM05935

## 产品信息

**Catalog Number:**  
CM05935

**CAS号:**  
870281-34-8


**分子式:**  
C<sub>21</sub>H<sub>16</sub>FN<sub>7</sub>O

**主要靶点:**  
PI3K

**主要通路:**  
PI3K/Akt/mTOR信号通路

**分子量:**  
401.4

**溶解度:**  
DMSO:73 mg/mL(181.9 mM)



## 靶点活性

p110 $\delta$ :14 nM

## 描述

Acalisib (GS-9820) is a potent and selective inhibitor of PI3K $\delta$ .

## 储存

Powder: -20°C for 3 years | In solvent: -80°C for 2 years

For technical support and original validation data for this product please contact

T: 027-87531629

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.