

Catalog Number: CM06692

产品信息

Catalog Number:
CM06692

CAS号:
1301761-96-5


分子式:
C₁₉H₂₃N₅O₄S

主要靶点:
Others

主要通路:
其他

分子量:
417.48

溶解度:
DMSO:5.5 mg/mL (13.17
mM),Need ultrasound



靶点活性

TNNI3K:25 nM

描述

GSK-114 is a selective inhibitor of TNNI3 interacting kinase (TNNI3K) with an IC₅₀ of 25 nM, 40-fold significant bias for TNNI3K compared to B-Raf with an IC₅₀ of 1000 nM, so GSK-114 has excellent Broad-spectrum kinase selectivity and sufficient oral exposure enable use in cellular and in vivo studies.

储存

Powder: -20°C for 3 years | In solvent: -80°C for 2 years