

For Research Use Only

BIM 23056 Acetate



www.ptgcn.com

Catalog Number: CM07397

产品信息

Catalog Number:
CM07397

CAS号:
T21944L

分子式:
 $C_{73}H_{85}N_{11}O_{11}$

主要靶点:
Others

主要通路:
其他

分子量:
1292.52

溶解度:

DMSO:60 mg/ml (46.42 mM)

CAS No.: T21944L

靶点活性

SSTR3:0.02 nM

描述

BIM 23056 Acetate is a potent and surmountable human recombinant sst5 receptor antagonist.

储存

Powder: -20°C for 3 years | In solvent: -80°C for 2 years

For technical support and original validation data for this product please contact

T: 027-87531629

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.