

Catalog Number: CM08263

### 产品信息

**Catalog Number:**  
CM08263

**CAS号:**  
300-84-5


**分子式:**  
C<sub>2</sub>H<sub>7</sub>NO<sub>2</sub>S

**主要靶点:**  
Chloride channel|Endogenous  
Metabolite

**主要通路:**  
代谢|离子通道

**分子量:**  
109.15

**溶解度:**  
H<sub>2</sub>O:150 mg/mL (1374.26 mM)



### 描述

Hypotaurine, an intermediate in taurine biosynthesis from cysteine in astrocytes, is an endogenous inhibitory amino acid of the glycine receptor.

### 储存

Powder: -20°C for 3 years | In solvent: -80°C for 2 years