

For Research Use Only

Enoxaparin sodium



www.ptgcn.com

Catalog Number: CM08572

产品信息

Catalog Number:
CM08572

CAS号:
679809-58-6

分子式:
 $C_{42}H_{59}N_3Na_4O_{35}S_2$

主要靶点:
Thrombin

主要通路:
蛋白酶体

分子量:
1322

溶解度:
DMSO:Insoluble

CAS No.: 679809-58-6

描述

Enoxaparin sodium is a low-molecular-weight heparin (LMWH) that has recently been approved by the US FDA for use in patients with medically managed ST-segment myocardial infarction (STEMI), or STEMI with subsequent percutaneous coronary intervention (PCI). It binds to and potentiates the action of antithrombin, and inhibits coagulation factors XIa, IXa, Xa and IIa (thrombin), thereby preventing formation of blood clots.

储存

Powder: -20°C for 3 years | In solvent: -80°C for 2 years

For technical support and original validation data for this product please contact

T: 027-87531629

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.