

For Research Use Only

Fas C-Terminal Tripeptide Acetate



www.ptgcn.com

Catalog Number: CM09734

产品信息

Catalog Number:
CM09734

CAS号:
T21724L


分子式:
 $C_{18}H_{33}N_3O_8$

主要靶点:
Fatty Acid Synthase|Others

主要通路:
代谢|其他

分子量:
419.47

溶解度:
DMSO:65 mg/ml (154.96 mM)



描述

Fas C-Terminal Tripeptide Acetate shows inhibitory activity of Fas/FAP-1 binding.

储存

Powder: -20°C for 3 years | In solvent: -80°C for 2 years

For technical support and original validation data for this product please contact

T: 027-87531629

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.