

For Research Use Only

Bimatoprost acid



www.ptgcn.com

Catalog Number: CM10319

产品信息

Catalog Number:
CM10319

CAS号:
38344-08-0


分子式:
C₂₃H₃₂O₅

主要靶点:
GPR

主要通路:
内分泌与激素|G蛋白偶联受体

分子量:
388.5

溶解度:
DMSO:65 mg/ml (167.31 mM)



描述

Bimatoprost acid is a metabolically stable analog of PGF₂ α and has a potential antagonistic activity for the FP receptor. It has a relative potency of 756% compared to PGF₂ α for binding to the FP receptor on ovine luteal cells.

储存

Powder: -20°C for 3 years | In solvent: -80°C for 2 years

For technical support and original validation data for this product please contact

T: 027-87531629

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.