

For Research Use Only

IL-4RA Protein, Mouse, Recombinant (hFc)



Catalog Number: CM10453

产品信息

Catalog Number:
CM10453

CAS号:
TMPJ-00387

分子式:
主要靶点:

分子量:
65-80 KDa, reducing conditions

CAS No.: TMPJ-00387

描述

Recombinant Mouse Interleukin-4 Receptor Subunit Alpha is produced by our Mammalian expression system and the target gene encoding Ile26-Arg233 is expressed with a human IgG1 Fc tag at the C-terminus.

储存

For technical support and original validation data for this product please contact

T: 027-87531629

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.