

For Research Use Only

Isoxsuprine hydrochloride



www.ptgcn.com

Catalog Number: CM11191

产品信息

Catalog Number:
CM11191

CAS号:
579-56-6


分子式:
C₁₈H₂₄ClNO₃

主要靶点:
Adrenergic Receptor|βGluR

主要通路:
G蛋白偶联受体|离子通道|神经科学

分子量:
337.84

溶解度:
DMSO:201.3 mM,H₂O:11.8 mM



靶点活性

myometrial beta-adrenergic receptor:13.65 μ M (Ki)|placental beta-adrenergic receptor:3.48 μ M (Ki)

描述

Isoxsuprine hydrochloride is a beta-adrenergic agonist that causes direct relaxation of uterine and vascular smooth muscle.

储存

Powder: -20°C for 3 years | In solvent: -80°C for 2 years

For technical support and original validation data for this product please contact

T: 027-87531629

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.