

For Research Use Only

# Isosorbide dinitrate



www.ptgcn.com

Catalog Number: CM12995

## 产品信息

**Catalog Number:**  
CM12995

**CAS号:**  
87-33-2


**分子式:**  
C<sub>6</sub>H<sub>8</sub>N<sub>2</sub>O<sub>8</sub>

**主要靶点:**  
Others|NO Synthase

**主要通路:**  
免疫与炎症|其他

**分子量:**  
236.14

**溶解度:**  
DMSO:199 mM



## 描述

Isosorbide dinitrate is a vasodilator used in the treatment of angina pectoris.

## 储存

Powder: -20°C for 3 years | In solvent: -80°C for 2 years

For technical support and original validation data for this product please contact

T: 027-87531629

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.