

For Research Use Only

C-MYC PEPTIDE EPI TOPE TFA



www.ptgcn.com

Catalog Number: CM13520

产品信息

Catalog Number:
CM13520

分子量:
1317.32

CAS号:
TP2311

分子式:
 $C_{53}H_{87}F_3N_{12}O_{23}$

Glu-Gln-Lys-Leu-Ile-Ser-Glu-Glu-Asp-Leu

主要靶点:
Others

主要通路:
其他

描述

C-MYC PEPTIDE EPI TOPE TFA, the product of the c-myc proto-oncogene, is a helix-loop-helix leucine zipper phosphoprotein that regulates gene transcription in cell proliferation, cell differentiation and apoptosis. This peptide is a human c-myc epitope.

储存

Powder: -20°C for 3 years | In solvent: -80°C for 2 years

For technical support and original validation data for this product please contact

T: 027-87531629

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.