

For Research Use Only  
Acetyltrialanine

Catalog Number: CM13625

产品信息

Catalog Number:  
CM13625


CAS号:  
19245-85-3

分子式:  
C<sub>11</sub>H<sub>19</sub>N<sub>3</sub>O<sub>5</sub>

主要靶点:

分子量:  
273.29

溶解度:  
DMSO:Soluble



描述

Acetyltrialanine is a dipeptide compound that binds at two sites on the Tb+3-pancreatic elastase complex and can be used as a nitrogen source.

储存

Powder: -20°C for 3 years | In solvent: -80°C for 2 years