

For Research Use Only

Fas Ligand Protein, Human, Recombinant (His)



Catalog Number: CM13893

产品信息

Catalog Number:
CM13893

分子量:
Approxiamtely 19.3 kDa

CAS No.:

CAS号:
分子式:
主要靶点:

描述

Fas Ligand, also known as FASLG and CD95L, is the ligand for FAS. It is a transmembrane protein which binds to TNFRSF6/FAS. Interaction of FAS with fas Ligand is critical in triggering apoptosis of some types of cells such as lymphocytes. Fas Ligand may be involved in cytotoxic T-cell mediated apoptosis and in T-cell development. TNFRSF6/FAS-mediated apoptosis may have a role in the induction of peripheral tolerance, in the antigen-stimulated suicide of mature T-cells, or both.

储存

For technical support and original validation data for this product please contact

T: 027-87531629

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.