

For Research Use Only

Rac1 Inhibitor W56 acetate(1095179-01-3 free base)



Catalog Number: CM14730

产品信息

Catalog Number:
CM14730

CAS号:
TP2131L

分子式:
 $C_{76}H_{121}N_{19}O_{25}S$

主要靶点:
Rho

主要通路:
细胞周期

分子量:
1732.95

溶解度:
DMSO:5 mM

Met-Val-Asp-Gly-Lys-Pro-Val-Asn-
Leu-Gly-Leu-Trp-Asp-Thr-Ala-Gly(acetate salt)

描述

Peptide comprising residues 45-60 of the guanine nucleotide exchange factor (GEF) recognition/activation site of Rac1; selectively inhibits Rac1 interaction with Rac1-specific GEFs TrioN, GEF-H1 and Tiam1.

储存

Powder: -20°C for 3 years | In solvent: -80°C for 2 years

For technical support and original validation data for this product please contact

T: 027-87531629

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.