

For Research Use Only

CEF20 acetate(153045-21-7 free base)



www.ptgcn.com

Catalog Number: CM15006

产品信息

Catalog Number:
CM15006

CAS号:
TP1687L

分子式:
 $C_{44}H_{78}N_{10}O_{14}S$

主要靶点:
Others

主要通路:
其他

分子量:
1003.22

溶解度:
DMSO:10 mM

H-Asn-Leu-Val-Pro-Met-Val-Ala-
Thr-Val-OH(acetate salt)

描述

CEF20 acetate is a CEF control peptide pool, is an HLA-A*0201-restricted epitope from cytomegalovirus pp65 (495-503).

储存

Powder: -20°C for 3 years | In solvent: -80°C for 2 years

For technical support and original validation data for this product please contact

T: 027-87531629

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.