

For Research Use Only

Boc-NH-PEG11-OH




www.ptgcn.com

Catalog Number: CM15329

产品信息

Catalog Number:
CM15329

分子量:
601.72



CAS号:
1556847-53-0

分子式:
 $C_{27}H_{55}NO_{13}$

主要靶点:
Others|PROTAC Linker

主要通路:
其他|PROTAC

体外活性

PROTACs contain two different ligands connected by a linker; one of them is a ligand for an E3 ubiquitin ligase and the other one is for the target protein. PROTACs exploit the intracellular ubiquitin-proteasome system to selectively degrade target proteins[1].

描述

Boc-NH-PEG11-OH is a PEG-based PROTAC linker that can be used in PROTAC synthesis.

储存

Powder: -20°C for 3 years | In solvent: -80°C for 2 years

For technical support and original validation data for this product please contact

T: 027-87531629

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.