

For Research Use Only

Disperse Fast Yellow 4K



www.ptgcn.com

Catalog Number: CM15792

产品信息

Catalog Number:
CM15792

CAS号:
6300-37-4


分子式:
C₁₉H₁₆N₄O

主要靶点:
Others

主要通路:
其他

分子量:
316.36

溶解度:
DMSO:3.16 mg/mL (10 mM),Need
ultrasound



描述

Disperse Fast Yellow 4K is a color additive used in dyes.

储存

Powder: -20°C for 3 years | In solvent: -80°C for 2 years

For technical support and original validation data for this product please contact

T: 027-87531629

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.