

For Research Use Only

# YAP-TEAD-IN-1 acetate



www.ptgcn.com

Catalog Number: CM15847

## 产品信息

**Catalog Number:**  
CM15847


**分子量:**  
2079.92

**CAS号:**  
TP2160L1

**分子式:**  
 $C_{95}H_{148}ClN_{23}O_{23}S_2$

**主要靶点:**  
YAP

**主要通路:**  
干细胞



## 体外活性

YAP-TEAD-IN-1 acetate (2  $\mu$ M) shows a higher binding affinity to TEAD1 (Kd = 15 nM) than YAP (50-171) (Kd = 40 nM)[1].

## 描述

YAP-TEAD-IN-1 acetate is an effective and competitive inhibitor of YAP-TEAD interaction with an IC50 of 25 nM.

## 储存

Powder: -20°C for 3 years | In solvent: -80°C for 2 years

For technical support and original validation data for this product please contact

T: 027-87531629

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.