

For Research Use Only

Mocravimod hydrochloride



www.ptgcn.com

Catalog Number: CM16909

产品信息

Catalog Number:
CM16909

CAS号:
509088-69-1


分子式:
C₂₄H₂₇Cl₂NO₃S

主要靶点:
LPL Receptor|Others

主要通路:
G蛋白偶联受体|其他

分子量:
480.45

溶解度:
DMSO:95 mg/mL (197.73 mM), Sonification is recommended.



体内活性

In mHC-disparate rat heart allografts, Mocravimod hydrochloride(orally; 0.1 and 1 mg/kg/day; for 100 days) prolongs graft survival and attenuates chronic rejection[1]. Mocravimod hydrochloride ameliorates atherosclerosis in LDL-R(-/-) Mice[2].

描述

Mocravimod hydrochloride is an effective and orally active agonist of sphingosine 1-phosphate receptor type 1(S1PR1).

储存

Powder: -20°C for 3 years | In solvent: -80°C for 2 years

For technical support and original validation data for this product please contact

T: 027-87531629

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.