

For Research Use Only

IL-6 Protein, Human, Recombinant



www.ptgcn.com

Catalog Number: CM18119

产品信息

Catalog Number:
CM18119

分子量:
Approxiamtely 20.95 kDa

CAS No.:

CAS号:
分子式:
主要靶点:

描述

Interleukin-6 (IL-6) is a multifunctional α -helical cytokine that regulates cell growth and differentiation of various tissues, which is known particularly for its role in the immune response and acute phase reactions. IL-6 protein is secreted by a variety of cell types including T cells and macrophages as a phosphorylated and variably glycosylated molecule. It exerts actions through its heterodimeric receptor composed of IL-6R that lacks the tyrosine/kinase domain and binds IL-6 with low affinity, and ubiquitously expressed glycoprotein 130 (gp130) that binds the IL-6. IL-6R complex with high affinity and thus transduces signals. IL-6 is also involved in hematopoiesis, bone metabolism, and cancer progression, and has been defined as an essential role in directing the transition from innate to acquired immunity. Cancer Immunotherapy/Immune Checkpoint Immunotherapy/Targeted Therapy

储存

For technical support and original validation data for this product please contact

T: 027-87531629

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.